



Glove Integrity Test System

Rechargeable Mobile Trolley

Rechargeable Mobile Trolley

Perfect Solution for Storing and Transporting Glove Integrity Testers

Tailin Rechargeable Mobile Trolley is a specialized device for use with Glove Integrity Testers. It makes it convenient to store and transport Glove Integrity Testers to the location where testing is required. Designed to make your work easier and more efficient.



Technical Parameters

Product	Tailin Rechargeable Mobile Trolley
Supporting Equipment	With 6 glove integrity testers
Charging Times	Provide 10 full charges
Human-Computer Interaction	Touch screen + 10.5" tablet
Transmission	WIFI
Dimension(mm)	1210×450×1085
Weight	≈85kg
Power Supply	220V, 50Hz
Power	0.15kW